

## Introduction and Background

Ten per cent of the worldwide adult population have some form of kidney disease, and every year millions die prematurely of complications related to Chronic Kidney Disease (CKD),<sup>1</sup> with 2.8 million End Stage Renal Disease (ESRD) patients routinely receiving dialysis.<sup>2</sup> The vast majority of patients dialyse 3 times weekly in-centre, with limited treatment alternatives. Quanta Dialysis Technologies has developed a high-performance, compact, easy-to-use haemodialysis (HD) system called SC+. The SC+ System is designed to bring flexibility to the dialysis community by allowing a greater choice in treatment frequency, setting location and regimen for dialysis patients. Quanta is conducting a Post Market Clinical Follow-up (PMCF) study to assess the safety, efficacy and usability of SC+. This analysis reports the results of 915 treatments completed between 30 June 2015 and 5 May 2017, and provides an update on the results previously reported.<sup>3</sup>

## Study Design

A prospective, multi-centre, open-label, single-arm observational study using SC+ was conducted with adult ESRD subjects receiving regular haemodialysis and weighing at least 50 kg. Each clinic's Principal Investigator (PI) screened the patients on a volunteer basis. Each dialysis therapy was performed using bicarbonate-buffered dialysate at a flow rate (Qd) of 500 ml/min. Vital signs were recorded nominally every 30 minutes, subject symptoms were monitored during every treatment and any adverse events were recorded. Blood samples for clinical laboratory evaluations of chemistry and haematology were taken pre- and post-dialysis in accordance with local practice, and subjects' mass was recorded before and after each treatment. Dialysis adequacy of SC+ was assessed by StdKt/V.

## Results

A total of 55 subjects with a mean age of 61 years were treated with SC+ on treatment regimens ranging from two times weekly to seven times per week. Subjects were screened by the PI at the clinic in line with the protocol requirements and stated intended use of SC+, with subjects' dry weights ranging from 50 to 157 kg. Table 1 summarises the patient characteristics of the study population. A total of 975 dialysis treatments were completed according to protocol during the period, for which 915 evaluable patient records are available. Clinical parameters are reported in Table 2. The mean duration of each treatment was 240.9 minutes, with a mean blood flow rate setting of 369.9 ml/min and dialysate flow rate of 500 ml/min. The range of blood pump speeds programmed on SC+ for the study population is displayed in Figure 1, with a study minimum of 220 ml/min and a maximum of 450 ml/min.

Only one ADE (Adverse Device Effect) was reported, where the device's safety system correctly detected a fault towards the end of treatment and subsequently instructed the operator to immediately end the session. As washback was not possible, the small volume of blood within the extra corporeal circuit was lost. This event was not classified as serious and the safety system of SC+ acted as intended to protect the patient. A total of 383 AE (Adverse Event) observations were recorded. With the exception of the single event noted above, all AEs related only to minor complications common in normal haemodialysis. The most frequent of these was hypotension, with 250 observations in 147 treatment sessions (16%) of treatment sessions, which is well within the normal range as documented in the literature.<sup>4</sup>

The association between ultrafiltration (UF) target and net fluid removal (NFR) achieved was assessed and is displayed in Figure 2. NFR achieved was determined by measuring the difference in patient mass before and after treatment (by using clinical scales) and corrected for priming, washback volumes, saline bolus and total oral intake. UF target matched with NFR achieved, with both parameters averaging 2.0 kg. The fit line for the data is well within the NFR error specified by International Standards for haemodialysis.<sup>5</sup> It should be noted that measurements of NFR obtained in clinical practice using weighing scales are prone to experimental error due to the influence of a number of variables that are difficult to control sufficiently in a real-world setting, such as intradialytic food and fluid intake and changes in articles of clothing and personal items.

Number of Subjects	55
Number of Evaluable Treatments	915
Date Range	30 June 2015 – 05 May 2017
Patient Age Range (years)	22 to 85; Mean 61
Patient Gender	32 Males, 23 Females
Dry Weight Range (Kg)	50 to 157; Mean 82.2
Treatments per Week	2/wk (n=1) ; 3/wk (n=48); 4/wk (n=3); 5/wk (n=2); 7/wk (n=1)
Anti-coagulant Type	Heparin (n=1); LMWH (n=44); Saline flush (n=1); None (n=9)
LMWH Total Dose (mg)	Dalteparin 0.45 – 5 (n = 9); Enoxaparin 2.5 – 60 (n = 28); Tinzaparin 0.45 – 5 (n = 7)
Needle Size (G)	14 Sharp (n=6); 14 Buttonhole (n=2); 15 Sharp (n=28); 15 Buttonhole (n=9); 16 Sharp (n=10)
Vascular Access	Radial AVF (n=14); Brachial Cephalic AVF (n =34); Brachial Cephalic AVG (n=3); Femoral Graft (n=4)
Blood Pump Speed (Qb) Range (ml/min)	220 – 450
Dialyser Size Range	FX-60 (n=200); FX-80 (n=189); FX-100 (n=272); Polyflux 170H (n=106); Polyflux 210H (n=148)
Dialysate Flow Rate (Qd) (ml/min)	500

Table 1: Baseline patient characteristics and dialysis treatment parameters for evaluable patients




Figure 1: Distribution of mean blood pump speed




Figure 2: Target ultrafiltration on SC+ plotted against the actual fluid removed by weight change pre- and post-treatment

## Clearance (StdKt/V)

Blood samples were collected monthly for each patient in accordance with individual clinic routines. Pre- and post-treatment blood analyses were available for 105 treatments within the study period and single pool Kt/V (spKt/V) for urea was calculated. These values were extrapolated to weekly standard urea Kt/V.<sup>6</sup> Figure 3 plots the StdKt/V for all available blood analyses. Subjects on a 5-times or 7-times weekly dialysis regimen displayed higher weekly clearance values, with a study maximum of 4.74 in a patient dialysing 7-times per week at home. The majority of data points were above the StdKt/V 2.1 weekly adequacy target, as recommended by Gotch.<sup>7</sup> Eleven data points did not reach this threshold. In 10 of these treatments the patient had significant residual renal function (RRF) and the dialysis prescriptions were adjusted accordingly. Notably, the patient on a two-times per week regimen (three hour treatment duration) had urine output of 1,000 ml/24h resulting in this low-dose treatment regimen. The remaining data point that did not reach the adequacy threshold was from a patient who had been unwell prior to dialysis involving hospital admission meaning that pre-treatment urea may have been low due to reduced food intake. Mean StdKt/V for patients without RRF was 2.47. Clinical performance was further demonstrated by a mean urea reduction ratio (URR) of 75.87% for the available blood data from patients without RRF, surpassing the recommended threshold of 65%.<sup>8</sup>




Figure 3: StdKt/V plotted for study population relative to the target of 2.1

## Technical Performance

The SC+ system consistently delivered key treatment parameters in line with the settings input by the operator, as shown in Table 2. In addition, deviation of dialysate temperature and conductivity from the target set point was negligible. Mean deviation of dialysate temperature was 0.1°C and mean deviation in conductivity was 0.01 mS/cm, demonstrating that SC+ maintains tight control over these parameters.

	Clinical Settings Mean ± SD	Clinical Measurements Mean ± SD
Mean Treatment Time (min)	240.9 ± 49.7 (target)	241.9 ± 43.0
Blood Pump speed (ml/min)	369.9 ± 46.8	370.0 ± 47.4
Dialysate Temperature (°C)	36.4 ± 0.4	36.4 ± 0.4
Dialysate Conductivity (mS/cm)	13.8 ± 0.1	13.7 ± 0.1
Target UF (kg)	2.0 ± 0.8	2.0 ± 0.8
Arterial Pressure (mmHg)	N/A	-127.8 ± 37.9
Venous Pressure (mmHg)	N/A	173.1 ± 32.4

Table 2: Clinical settings and machine measurements for evaluable treatments conducted

## Discussion and Conclusion

This PMCF study supports the safety and performance of the SC+ Haemodialysis System with acceptable clinical and technical performance in a diverse sample of haemodialysis subjects. There was only one ADE, which was minor, with the safety system of SC+ acting as intended to protect the patient from harm. All other AEs were typical of dialysis treatments and not device-related. Clinical performance was demonstrated by subjects consistently attaining a weekly StdKt/V of greater than 2.1 as recommended by Gotch as a weekly dialysis adequacy.<sup>7</sup> Clinical performance was further demonstrated by a mean study URR of 74.5%, surpassing the recommended guidelines of 65% for 3-times weekly dialysis. UF performance was acceptable, with actual measured UF closely matching the target UF. SC+ performed consistently across all treatments, with results overall demonstrating that patients treated with SC+ achieve adequate clearance using standard treatment regimens traditionally used in haemodialysis. Results further show that SC+ can be used safely and effectively in a population of adult haemodialysis patients representative of those found in a typical dialysis clinic, with no significant safety issues identified.

1. Couser WG, Remuzzi G, Mendis S, Tonelli M. The contribution of chronic kidney disease to the global burden of major noncommunicable diseases. *Kidney Int*. Dec 2011; 80(12): 1258-70.

2. Fresenius Medical Care Annual Report 2015

3. Poster presented at the Annual Dialysis Conference, Long Beach, California, March 2017

4. Levy J. *Oxford Handbook of Dialysis*, Third Edition

5. IEC60601-2-16:2012, available at <https://webstore.iec.ch/publication/2619>, accessed 3 March 2017

6. Daugirdas JT. Second generation logarithmic estimates of single-pool variable volume Kt/V: an analysis of error. *J Am Soc Nephrol* 1993; 4: 1205-13

7. Gotch FA. The current place of urea kinetic modelling with respect to different dialysis modalities. *Nephrol Dial Transplant* 1998; 13 (Suppl 6):S10-4

8. UK Renal Association Guidelines: <http://www.renal.org/docs/default-source/guidelines-resources/old-guidelines/haemodialysis-5th-editionccc0fa231181561659443ff00014d4d8.pdf?sfvrsn=2>, accessed on 23 May 2017